

# Dispenser Selection Chart

Reagent	Dispensette® S	Dispensette® Organic	Reagent	Dispensette® S	Dispensette® Organic	Reagent	Dispensette® S	Dispensette® Organic
Acetaldehyde	+	+	Cyclohexanone	+	+	Methylene chloride		+
Acetic acid (glacial), 100%	+	+	Cyclopentane		+	Mineral oil (Engine oil)	+	+
Acetic acid, ≤ 96%	+	+	Decane	+	+	Monochloroacetic acid	+	+
Acetic anhydride		+	1-Decanol	+	+	Nitric acid, ≤ 30%	+	+
Acetone	+	+	Dibenzyl ether	+	+	Nitric acid, 30-70% */ **		+
Acetonitrile	+	+	Dichloroacetic acid		+	Nitrobenzene	+	+
Acetophenone		+	Dichlorobenzene	+	+	Oleic acid	+	+
Acetyl chloride		+	Dichloroethane		+	Oxalic acid	+	
Acetylacetone	+	+	Dichloroethylene		+	n-Pentane		+
Adipic acid	+		Dichloromethane		+	Peracetic acid		+
Allyl alcohol	+	+	Diesel oil (Heating oil), bp 250-350 °C		+	Perchloric acid	+	+
Aluminium chloride	+		Diethanolamine	+	+	Perchloroethylene		+
Amino acids	+		Diethyl ether		+	Petroleum, bp 180-220 °C		+
Ammonia, ≤ 20%	+	+	Diethylamine	+	+	Petroleum ether, bp 40-70 °C		+
Ammonia, 20-30%		+	1,2 Diethylbenzene	+	+	Phenol	+	+
Ammonium chloride	+		Diethylene glycol	+	+	Phenylethanol	+	+
Ammonium fluoride	+		Dimethyl sulfoxide (DMSO)	+	+	Phenylhydrazine	+	+
Ammonium sulfate	+		Dimethylaniline	+		Phosphoric acid, ≤ 85%	+	+
n-Amyl acetate	+	+	Dimethylformamide (DMF)	+	+	Phosphoric acid, 85% + Sulfuric acid, 98%, 1:1	+	+
Amyl alcohol (Pentanol)	+	+	1,4 Dioxane		+	Piperidine	+	+
Amyl chloride (Chloropentane)		+	Diphenyl ether	+	+	Potassium chloride	+	
Aniline	+	+	Essential oil		+	Potassium dichromate	+	
Barium chloride	+		Ethanol	+	+	Potassium hydroxide	+	
Benzaldehyde	+	+	Ethanolamine	+	+	Potassium permanganate	+	
Benzene (Benzol)	+	+	Ethyl acetate	+	+	Propionic acid	+	+
Benzene (Petroleum benzin), bp 70-180 °C		+	Ethylbenzene		+	Propylene glycol (Propanediol)	+	+
Benzoyl chloride	+	+	Ethylene chloride		+	Pyridine	+	+
Benzyl alcohol	+	+	Fluoroacetic acid		+	Pyruvic acid	+	+
Benzylamine	+	+	Formaldehyde, ≤ 40%	+		Salicylaldehyde	+	+
Benzylchloride	+	+	Formamide	+	+	Scintillation fluid	+	+
Boric acid, ≤ 10%	+	+	Formic acid, ≤ 100%		+	Silver acetate	+	
Bromobenzene	+	+	Glycerol	+	+	Silver nitrate	+	
Bromonaphthalene	+	+	Glycol (Ethylene glycol)	+	+	Sodium acetate	+	
Butanediol	+	+	Glycolic acid, ≤ 50%	+		Sodium chloride	+	
1-Butanol	+	+	Heating oil (Diesel oil), bp 250-350 °C		+	Sodium dichromate	+	
n-Butyl acetate	+	+	Heptane		+	Sodium fluoride	+	
Butyl methyl ether	+	+	Hexane		+	Sodium hydroxide, ≤ 30%	+	
Butylamine	+	+	Hexanoic acid	+	+	Sodium hypochlorite	+	
Butyric acid	+	+	Hexanol	+	+	Sulfuric acid, ≤ 98%	+	+
Calcium carbonate	+		Hydriodic acid, ≤ 57% **	+	+	Tartaric acid	+	
Calcium chloride	+		Hydrobromic acid	+	+	Tetrachloroethylene		+
Calcium hydroxide	+		Hydrochloric acid, ≤ 20%	+	+	Tetrahydrofuran (THF) */ **		+
Calcium hypochlorite	+		Hydrochloric acid, 20-37% **		+	Tetramethylammonium hydroxide	+	
Carbon tetrachloride		+	Hydrogen peroxide, ≤ 35%		+	Toluene		+
Chloro naphthalene	+	+	Isoamyl alcohol	+	+	Trichloroacetic acid		+
Chloroacetaldehyde, ≤ 45%	+	+	Isobutanol	+	+	Trichlorobenzene		+
Chloroacetic acid	+	+	Isooctane		+	Trichloroethane		+
Chloroacetone	+	+	Isopropanol (2-Propanol)	+	+	Trichloroethylene		+
Chlorobenzene	+	+	Isopropyl ether	+	+	Trichlorotrifluoro ethane		+
Chlorobutane	+	+	Lactic acid	+		Triethanolamine	+	+
Chloroform		+	Methanol	+	+	Triethylene glycol	+	+
Chlorosulfonic acid		+	Methoxybenzene	+	+	Trifluoro ethane		+
Chromic acid, ≤ 50%	+	+	Methyl benzoate	+	+	Trifluoroacetic acid (TFA)		+
Chromosulfuric acid	+		Methyl butyl ether	+	+	Turpentine		+
Copper sulfate	+		Methyl ethyl ketone	+	+	Urea	+	
Cresol		+	Methyl formate	+	+	Xylene		+
Cumene (Isopropyl benzene)	+	+	Methyl propyl ketone	+	+	Zinc chloride, ≤ 10%	+	
Cyclohexane		+				Zinc sulfate, ≤ 10%	+	

The above recommendations reflect testing completed prior to publication. Always follow instructions in the operating manual of the instrument as well as the reagent manufacturer's specifications. In addition to these chemicals, a variety of organic and inorganic saline solutions (e.g., biological buffers), biological detergents and media for cell culture can be dispensed. Should you require information on chemicals not listed, please feel free to contact BRAND. Status as of: 01/24/14

\* use ETFE/PTFE bottle adapter

\*\* use PTFE seal for valve block

For dispensing HF, we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring. Please find further product information on [www.brand.de](http://www.brand.de)

