



Buy one  
get the  
second one  
absolutely  
**free\***

## Transferpette® S -8

Don't miss this offer for multichannel micropipettes  
– only valid until October 31st, 2020.

**BRAND. For lab. For life.®**

### Get to know the new Transferpette® S -8

The Transferpette® S air displacement pipette allows you to work efficiently and ergonomically with both small and large volumes. Low operating forces and a short stroke ensure relaxed pipetting, even during long pipetting series. Due to its perfect ergonomics, the Transferpette® S provides a comfortable grip in any hand position, whether you are right- or left-handed, or have big or small hands. Changing your grip is not necessary, since the microliter pipette provides true one-hand operation: set the volume, pipette, eject the tip.

### Your benefits:

- + Efficient: one-hand operation
- + Flexible: ergonomic design for all hand sizes
- + Reliable: volume-change protection prevents accidental adjustments
- + Free ultra-low retention tips with your purchase.

Learn more in our video:



# Buy one – get one free promotion Transferpette® S -8

Validity period: February 1st - October 31st, 2020

## Buy one - get one free\*!

Just fill in the form and send it to  
customersupport@brand.co.in



o Yes, I would like to participate and order the following quantities as per the BRAND Scientific Equipment Private Ltd. General Terms and Conditions.

### Micropipette, Transferpette® S -8

Capacity, µl (color-coded)	Description	Cat. No.	Quantity**	Price ₹*	Total price ₹*
0.5 – 10	Transferpette® S-8, M8-10	705900		67,953	
5 – 50	Transferpette® S-8, M8-50	705906		67,953	
10 – 100	Transferpette® S-8, M8-100	705908		67,953	
20 – 200	Transferpette® S-8, M8-200	705910		67,953	
30 – 300	Transferpette® S-8, M8-300	705912		67,953	

**Total:**

Date

Signature

#### Address:

.....  
First name

.....  
Last name

.....  
Institute/Company

.....  
Street

.....  
Postal code and city

.....  
E-Mail

.....  
Tel./Fax

\* Only valid during promotion period. Customer will receive one instrument free from any of the volume ranges mentioned against each pc ordered. Price mentioned are on FOR basis.

All additional taxes and duties namely, GST & other levies (if any) wherever applicable will be charged extra.

\*\* Quantity ordered - same quantity will be additionally delivered free of charge.

Transferpette®, BRAND®, BRAND. For lab. For life.®, as well as the figurative marks depicted here and the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners.

BRAND Scientific Equipment Pvt. Ltd. | 303, 3rd Floor, 'C' Wing | Delphi, Hiranandani Business Park

Powai Mumbai | 400 076 (India)

T +91 22 42957790 | F +91 22 42957791 | customerservice@brand.co.in | www.brand.co.in

