

Test our new sealing film Seal-R-film™ and order through your distributor.



Seal-R-film™

For easy and efficient sealing, masking, and insulating.



Sealing film Seal-R-film™

- Protects against evaporation, contamination, and leakage
- Also fits irregular shapes
- Simply cut from the roll and seal your vessel

BRAND. For lab. For life.




NOW!

Seal-R-film™ Try and Buy

The all-purpose solution for sealing in the laboratory

Beakers, flasks, and other vessels in the laboratory frequently have to be stored temporarily. BRAND Seal-R-film™ sealing film ensures that your samples remain protected against evaporation, contamination, and leakage. The stretchable film is easy to handle, as it only adheres once it has been stretched. Otherwise, it adheres neither to itself nor to gloves. The film's resistance to a large number of chemicals makes it fit for universal use in the laboratory.



Easy to use

Simply take Seal-R-film™ from the roll, remove the paper, and pull it over the area to be sealed while stretching the film uniformly. The optional cutter allows the clean storage and convenient cutting of the sealing film.

Request a sample and order Seal-R-film™ now
Simply place your order through your distributor.

Seal-R-film™

Description	Pack of	Cat. No.
Width 100 mm, Length 38 m	12 rolls	7016 20
Width 50 mm, Length 75 m	24 rolls	7016 22
Width 100 mm, Length 75 m	12 rolls	7016 24
Width 500 mm, Length 15 m	6 rolls	7016 26
Cutter	1 piece	7016 50



BRAND®, BRAND. For lab. For life.® as well as the logo design mark depicted are trademarks of BRAND GMBH + CO KG, Germany. Other brands and trademarks as referenced in the current printed matter are property of their respective third-party owners. California Residents: For more information concerning California Proposition 65, please refer to www.brand.de/calprop65. Subject to technical modification without notice. Errors excepted.

BRAND Scientific Equipment Pvt. Ltd. · 303, 3rd Floor, 'C' Wing, Delphi · Hiranandani Business Park, Powai
Mumbai-400076(India) · Tel: +91 2242957790 · Fax: +91 2242957791 · customersupport@brand.co.in · www.brand.co.in

BRAND GMBH + CO KG · P.O. Box 1155 · 97861 Wertheim · Germany
Tel.: +49 9342 808-0 · Fax: +49 9342 808-98000 · info@brand.de · www.brand.de

